



SH3RF3 Antibody

Product Code	CSB-PA009465
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEJ3
Immunogen	Synthesized peptide derived from the C-terminal region of Human SH3RF3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SH3RF3; POSH2; SH3MD4; SH3 domain-containing RING finger protein 3; Plenty of SH3s 2; SH3 multiple domains protein 4
Immunogen Species	Homo sapiens (Human)
Target Names	SH3RF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.