



# IFNAR1 Antibody

<b>Product Code</b>	CSB-PA009473
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17181
<b>Immunogen</b>	Synthesized peptide derived from Human IFN- $\alpha$ / $\beta$ R $\alpha$ around the non-phosphorylation site of Y466.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	IFNAR1; IFNAR; Interferon alpha/beta receptor 1; IFN-R-1; IFN-alpha/beta receptor 1; Cytokine receptor class-II member 1; Cytokine receptor family 2 member 1; CRF2-1; Type I interferon receptor 1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IFNAR1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.