




OR2A1; Antibody

Product Code	CSB-PA009480
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGT9
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 2A42.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR2A1; OR2A42; Olfactory receptor 2A1/2A42; Olfactory receptor OR7-16; Olfactory receptor OR7-19
Immunogen Species	Homo sapiens (Human)
Target Names	OR2A1

Image



Western Blot analysis of COLO cells using Olfactory receptor 2A42 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.