




# OR2A1; Antibody

<b>Product Code</b>	CSB-PA009480
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8NGT9
<b>Immunogen</b>	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 2A42.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	OR2A1; OR2A42; Olfactory receptor 2A1/2A42; Olfactory receptor OR7-16; Olfactory receptor OR7-19
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	OR2A1

## Image



Western Blot analysis of COLO cells using Olfactory receptor 2A42 Polyclonal Antibody