



Phospho-IGF2R (S2409) Antibody

Product Code	CSB-PA009491
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11717
Immunogen	Synthesized peptide derived from Human IGF-IIR around the phosphorylation site of S2409.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IGF2R; MPRI; Cation-independent mannose-6-phosphate receptor; CI Man-6-P receptor; CI-MPR; M6PR; 300 kDa mannose 6-phosphate receptor; MPR 300; Insulin-like growth factor 2 receptor; Insulin-like growth factor II receptor; IGF-II receptor;
Immunogen Species	Homo sapiens (Human)
Target Names	IGF2R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.