



Phospho-IGFBP3 (S183) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA009507 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P17936 |
| Immunogen | Synthesized peptide derived from Human IGFBP3 around the phosphorylation site of S183. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | IGFBP3; IBP3; Insulin-like growth factor-binding protein 3; IBP-3; IGF-binding protein 3; IGFBP-3 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IGFBP3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |