



IL2RA Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA009543 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P01589 |
| Immunogen | Synthesized peptide derived from Human IL-2R α around the non-phosphorylation site of S268. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | IL2RA; Interleukin-2 receptor subunit alpha; IL-2 receptor subunit alpha; IL-2-RA; IL-2R subunit alpha; IL2-RA; TAC antigen; p55; CD antigen CD25 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IL2RA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |