




Phospho-IL2RB (Y364) Antibody

Product Code	CSB-PA009547
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14784
Immunogen	Synthesized peptide derived from Human IL-2R β around the phosphorylation site of Y364.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL2RB; Interleukin-2 receptor subunit beta; IL-2 receptor subunit beta; IL-2R subunit beta; IL-2RB; High affinity IL-2 receptor subunit beta; p70-75; p75; CD antigen CD122
Immunogen Species	Homo sapiens (Human)
Target Names	IL2RB

Image



Western Blot analysis of 293T cells using Phospho-IL-2R β (Y364) Polyclonal Antibody