



IL4R Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA009556 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P24394 |
| Immunogen | Synthesized peptide derived from the Internal region of Human IL-4R α . |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IL4R |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |