




IL6ST Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA009560 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P40189 |
| Immunogen | Synthesized peptide derived from Human CD130 around the non-phosphorylation site of S782. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | IL6ST; Interleukin-6 receptor subunit beta; IL-6 receptor subunit beta; IL-6R subunit beta; IL-6R-beta; IL-6RB; CDw130; Interleukin-6 signal transducer; Membrane glycoprotein 130; gp130; Oncostatin-M receptor subunit alpha; CD antigen CD130 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IL6ST |

Image



Western Blot analysis of 293T COLO cells using CD130 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.