




Phospho-CXCR2 (S347) Antibody

Product Code	CSB-PA009566
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25025
Immunogen	Synthesized peptide derived from Human IL-8R β around the phosphorylation site of S347.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CXCR2; IL8RB; C-X-C chemokine receptor type 2; CXC-R2; CXCR-2; CDw128b; GRO/MGSA receptor; High affinity interleukin-8 receptor B; IL-8R B; IL-8 receptor type 2; CD antigen CD182
Immunogen Species	Homo sapiens (Human)
Target Names	CXCR2

Image



Western Blot analysis of 3T3 cells using Phospho-IL-8R β (S347) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.