




IL12A Antibody

Product Code	CSB-PA009581
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29459
Immunogen	Synthesized peptide derived from the Internal region of Human IL-12A.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ICC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ICC:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL12A; NKSF1; Interleukin-12 subunit alpha; IL-12A; Cytotoxic lymphocyte maturation factor 35 kDa subunit; CLMF p35; IL-12 subunit p35; NK cell stimulatory factor chain 1; NKSF1
Immunogen Species	Homo sapiens (Human)
Target Names	IL12A

Image



Western Blot analysis of HeLa cells using IL-12A Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.