



GPR153 Antibody

Product Code	CSB-PA009672
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NV75
Immunogen	Synthesized peptide derived from the Internal region of Human GPR153.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR153; PGR1; Probable G-protein coupled receptor 153; G-protein coupled receptor PGR1
Immunogen Species	Homo sapiens (Human)
Target Names	GPR153
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.