



OR8H3 Antibody

Product Code	CSB-PA009734
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N146
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 8H3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR8H3; Olfactory receptor 8H3; Olfactory receptor OR11-172
Immunogen Species	Homo sapiens (Human)
Target Names	OR8H3