



OR5M1/OR5M10 Antibody

Product Code	CSB-PA009740
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGP8/Q6IEU7
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 5M1/5M10.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR5M1; Olfactory receptor 5M1; OST050; Olfactory receptor OR11-208; OR5M10; Olfactory receptor 5M10; Olfactory receptor OR11-207
Immunogen Species	Homo sapiens (Human)
Target Names	OR5M1/OR5M10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.