




OR5B12 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA009742 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96R08 |
| Immunogen | Synthesized peptide derived from the Internal region of Human Olfactory receptor 5B12. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | OR5B12; OR5B12P; OR5B16; Olfactory receptor 5B12; Olfactory receptor 5B16; Olfactory receptor OR11-241 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | OR5B12 |

Image



Western Blot analysis of A549 cells using Olfactory receptor 5B12 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.