



OR13C2/OR13C9 Antibody

Product Code	CSB-PA009783
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGS9/Q8NGT0
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor O13C2/9.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR13C2; Olfactory receptor 13C2; Olfactory receptor OR9-12; OR13C9; Olfactory receptor 13C9; Olfactory receptor OR9-13
Immunogen Species	Homo sapiens (Human)
Target Names	OR13C2/OR13C9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.