




LIPC Antibody

Product Code	CSB-PA009836
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11150
Immunogen	Synthesized peptide derived from the Internal region of Human Hepatic Lipase.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LIPC; HTGL; Hepatic triacylglycerol lipase; HL; Hepatic lipase; Lipase member C
Immunogen Species	Homo sapiens (Human)
Target Names	LIPC

Image



Western Blot analysis of A549 cells using Hepatic Lipase Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.