




LIPE Antibody

Product Code	CSB-PA009841
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05469
Immunogen	Synthesized peptide derived from Human HSL around the non-phosphorylation site of S855.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LIPE; Hormone-sensitive lipase; HSL
Immunogen Species	Homo sapiens (Human)
Target Names	LIPE

Image



Western Blot analysis of HuvEc cells using HSL Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.