



MAP4 Antibody

Product Code	CSB-PA010006
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27816
Immunogen	Synthesized peptide derived from Human MAP-4 around the non-phosphorylation site of S696.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP4; Microtubule-associated protein 4; MAP-4
Immunogen Species	Homo sapiens (Human)
Target Names	MAP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.