




Phospho-TFRC (S24) Antibody

Product Code	CSB-PA010092
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02786
Immunogen	Synthesized peptide derived from Human CD71 around the phosphorylation site of S24.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD antigen CD71
Immunogen Species	Homo sapiens (Human)
Target Names	TFRC

Image



Western Blot analysis of 293 cells using Phospho-CD71 (S24) Polyclonal Antibody