



Phospho-MCM4 (S54) Antibody

Product Code	CSB-PA010108
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33991
Immunogen	Synthesized peptide derived from Human MCM4 around the phosphorylation site of S54.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MCM4; CDC21; DNA replication licensing factor MCM4; CDC21 homolog; P1-CDC21
Immunogen Species	Homo sapiens (Human)
Target Names	MCM4