



Phospho-MDM2 (S186/S188) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA010119 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q00987 |
| Immunogen | Synthesized peptide derived from Human MDM2 around the phosphorylation site of S186/S188. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | IF, ELISA; IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | MDM2; E3 ubiquitin-protein ligase Mdm2; Double minute 2 protein; Hdm2; Oncoprotein Mdm2; p53-binding protein Mdm2 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MDM2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |