




MDM4 Antibody

Product Code	CSB-PA010135
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15151
Immunogen	Synthesized peptide derived from Human MDMX around the non-phosphorylation site of S367.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MDM4; MDMX; Protein Mdm4; Double minute 4 protein; Mdm2-like p53-binding protein; Protein Mdmx; p53-binding protein Mdm4
Immunogen Species	Homo sapiens (Human)
Target Names	MDM4

Image



Western Blot analysis of 293T cells using MDMX Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.