



Phospho-MAP3K9/MAP3K10 (T312/266) Antibody

Product Code	CSB-PA010197
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80192/Q02779
Immunogen	Synthesized peptide derived from Human MLK1/2 around the phosphorylation site of T312/266.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP3K9; MLK1; PRKE1; Mitogen-activated protein kinase kinase kinase 9; Mixed lineage kinase 1; MAP3K10; MLK2; MST; Mitogen-activated protein kinase kinase kinase 10; Mixed lineage kinase 2; Protein kinase MST
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K9/MAP3K10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.