




ACHE Antibody

Product Code	CSB-PA010202
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22303
Immunogen	Synthesized peptide derived from the Internal region of Human AChE.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACHE; Acetylcholinesterase; AChE
Immunogen Species	Homo sapiens (Human)
Target Names	ACHE

Image



Western Blot analysis of MCF7 cells using AChE Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.