




MNDA Antibody

Product Code	CSB-PA010214
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41218
Immunogen	Synthesized peptide derived from the C-terminal region of Human MNDA.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MNDA; Myeloid cell nuclear differentiation antigen
Immunogen Species	Homo sapiens (Human)
Target Names	MNDA

Image



Western Blot analysis of 3T3 cells using MNDA Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.