



Phospho-EIF2AK4 (T899) Antibody

Product Code	CSB-PA010224
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2K8
Immunogen	Synthesized peptide derived from Human GCN2 around the phosphorylation site of T899.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EIF2AK4; GCN2; KIAA1338; Eukaryotic translation initiation factor 2-alpha kinase 4; GCN2-like protein
Immunogen Species	Homo sapiens (Human)
Target Names	EIF2AK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.