




OR13G1 Antibody

Product Code	CSB-PA010229
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGZ3
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 13G1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR13G1; Olfactory receptor 13G1; Olfactory receptor OR1-37
Immunogen Species	Homo sapiens (Human)
Target Names	OR13G1

Image



Western Blot analysis of 293 cells using Olfactory receptor 13G1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.