



OR14J1 Antibody

Product Code	CSB-PA010248
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGF5
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 5U1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR14J1; OR5U1; Olfactory receptor 14J1; Hs6M1-28; Olfactory receptor 5U1; Olfactory receptor OR6-25
Immunogen Species	Homo sapiens (Human)
Target Names	OR14J1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.