




MSN/RDX/EZR Antibody

Product Code	CSB-PA010258
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26038/P35241/P15311
Immunogen	Synthesized peptide derived from Human Moesin/Ezrin/Radixin around the non-phosphorylation site of T558.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MSN; Moesin; Membrane-organizing extension spike protein; RDX; Radixin; EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81
Immunogen Species	Homo sapiens (Human)
Target Names	MSN/RDX/EZR

Image



Western Blot analysis of 3T3 cells using Moesin/Ezrin/Radixin Polyclonal Antibody