




Phospho-MST1R (S1394) Antibody

Product Code	CSB-PA010260
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04912
Immunogen	Synthesized peptide derived from Human Ron around the phosphorylation site of S1394.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MST1R; PTK8; RON; Macrophage-stimulating protein receptor; MSP receptor; CDw136; Protein-tyrosine kinase 8; p185-Ron; CD antigen CD136
Immunogen Species	Homo sapiens (Human)
Target Names	MST1R

Image



Western Blot analysis of 3T3 cells using Phospho-Ron (S1394) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.