



MX2 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA010290 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P20592 |
| Immunogen | Synthesized peptide derived from the Internal region of Human MxB. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA; ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | MX2; Interferon-induced GTP-binding protein Mx2; Interferon-regulated resistance GTP-binding protein MxB; Myxovirus resistance protein 2; p78-related protein |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MX2 |