



MX2 Antibody

Product Code	CSB-PA010290
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20592
Immunogen	Synthesized peptide derived from the Internal region of Human MxB.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA; ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MX2; Interferon-induced GTP-binding protein Mx2; Interferon-regulated resistance GTP-binding protein MxB; Myxovirus resistance protein 2; p78-related protein
Immunogen Species	Homo sapiens (Human)
Target Names	MX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.