



ZFH3 Antibody

Product Code	CSB-PA020003
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15911
Immunogen	Synthesized peptide derived from the Internal region of Human ATBF1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZFH3; ATBF1; Zinc finger homeobox protein 3; AT motif-binding factor; AT-binding transcription factor 1; Alpha-fetoprotein enhancer-binding protein; Zinc finger homeodomain protein 3; ZFH-3
Immunogen Species	Homo sapiens (Human)
Target Names	ZFH3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.