



Phospho-MYLK (Y464) Antibody

Product Code	CSB-PA020005
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15746
Immunogen	Synthesized peptide derived from Human MYLK around the phosphorylation site of Y464.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MYLK; MLCK; MLCK1; MYLK1; Myosin light chain kinase; smooth muscle; MLCK; smMLCK; Kinase-related protein; KRP; Telokin
Immunogen Species	Homo sapiens (Human)
Target Names	MYLK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.