



MYO5A Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA020010 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y4I1 |
| Immunogen | Synthesized peptide derived from the C-terminal region of Human Myosin VA. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | MYO5A; MYH12; Unconventional myosin-Va; Dilute myosin heavy chain; non-muscle; Myosin heavy chain 12; Myosin-12; Myoxin |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MYO5A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |