



NDUFC2 Antibody

Product Code	CSB-PA020040
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95298
Immunogen	Synthesized peptide derived from the Internal region of Human NDUFC2.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NDUFC2; HLC1; NADH dehydrogenase [ubiquinone] 1 subunit C2; Complex I-B14.5b; CI-B14.5b; Human lung cancer oncogene 1 protein; HLC-1; NADH-ubiquinone oxidoreductase subunit B14.5b
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.