




ATM Antibody

Product Code	CSB-PA020042
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13315
Immunogen	Synthesized peptide derived from Human Atm around the non-phosphorylation site of S1981.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ATM; Serine-protein kinase ATM; Ataxia telangiectasia mutated; A-T mutated
Immunogen Species	Homo sapiens (Human)
Target Names	ATM

Image



Western Blot analysis of MOUSE-BRAIN K562 cells using Atm Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.