



NDUFS6 Antibody

Product Code	CSB-PA020046
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75380
Immunogen	Synthesized peptide derived from the C-terminal region of Human NDUFS6.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NDUFS6; NADH dehydrogenase [ubiquinone] iron-sulfur protein 6; mitochondrial; Complex I-13kD-A; CI-13kD-A; NADH-ubiquinone oxidoreductase 13 kDa-A subunit
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFS6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.