



NDUFV3 Antibody

Product Code	CSB-PA020047
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56181
Immunogen	Synthesized peptide derived from the Internal region of Human NDUFV3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NDUFV3; NADH dehydrogenase [ubiquinone] flavoprotein 3; mitochondrial; Complex I-9kD; CI-9kD; NADH-ubiquinone oxidoreductase 9 kDa subunit; Renal carcinoma antigen NY-REN-4
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFV3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.