




Phospho-NF2 (S518) Antibody

Product Code	CSB-PA020074
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35240
Immunogen	Synthesized peptide derived from Human NF2 around the phosphorylation site of S518.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NF2; SCH; Merlin; Moesin-ezrin-radixin-like protein; Neurofibromin-2; Schwannomerlin; Schwannomin
Immunogen Species	Homo sapiens (Human)
Target Names	NF2

Image



Western blot analysis of HELA using p-NF2 (S518) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.