



Phospho-NFATC1 (S294) Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA020078 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O95644 |
| Immunogen | Synthesized peptide derived from Human NFATc1 around the phosphorylation site of S294. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | NFATC1; NFAT2; NFATC; Nuclear factor of activated T-cells; cytoplasmic 1; NF-ATc1; NFATc1; NFAT transcription complex cytosolic component; NF-ATc; NFATc |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | NFATC1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |