



Phospho-NFKBIA (Y305) Antibody

Product Code	CSB-PA020112
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25963
Immunogen	Synthesized peptide derived from Human IκB-α around the phosphorylation site of Y305.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NFKBIA; IKBA; MAD3; NFKBI; NF-kappa-B inhibitor alpha; I-kappa-B-alpha; IκB-alpha; IkappaBalphα; Major histocompatibility complex enhancer-binding protein MAD3
Immunogen Species	Homo sapiens (Human)
Target Names	NFKBIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.