



# NFKBIA Antibody

<b>Product Code</b>	CSB-PA020113
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P25963
<b>Immunogen</b>	Synthesized peptide derived from Human IκB-α around the non-phosphorylation site of Y305.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	NFKBIA; IKBA; MAD3; NFKBI; NF-kappa-B inhibitor alpha; I-kappa-B-alpha; IκB-alpha; IkappaBalpa; Major histocompatibility complex enhancer-binding protein MAD3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	NFKBIA
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.