




Phospho-NFKBIB (T19) Antibody

Product Code	CSB-PA020117
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15653
Immunogen	Synthesized peptide derived from Human IκB-β around the phosphorylation site of T19.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NFKBIB; IKBB; TRIP9; NF-kappa-B inhibitor beta; NF-kappa-BIB; I-kappa-B-beta; IκB-B; IκB-beta; IκappaBbeta; Thyroid receptor-interacting protein 9; TR-interacting protein 9; TRIP-9
Immunogen Species	Homo sapiens (Human)
Target Names	NFKBIB

Image



Western Blot analysis of 293 cells using Phospho-IκB-β (T19) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.