




Phospho-ACO1 (S711) Antibody

Product Code	CSB-PA020121
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21399
Immunogen	Synthesized peptide derived from Human IRP-1 around the phosphorylation site of S711.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACO1; IREB1; Cytoplasmic aconitate hydratase; Aconitase; Citrate hydro-lyase; Ferritin repressor protein; Iron regulatory protein 1; IRP1; Iron-responsive element-binding protein 1; IRE-BP 1
Immunogen Species	Homo sapiens (Human)
Target Names	ACO1

Image



Western Blot analysis of 293 cells using Phospho-IRP-1 (S711) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.