



NKTR Antibody

Product Code	CSB-PA020126
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30414
Immunogen	Synthesized peptide derived from the Internal region of Human NK-TR.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NKTR; NK-tumor recognition protein; NK-TR protein; Natural-killer cells cyclophilin-related protein
Immunogen Species	Homo sapiens (Human)
Target Names	NKTR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.