





# NOS3 Antibody

<b>Product Code</b>	CSB-PA020137
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P29474
<b>Immunogen</b>	Synthesized peptide derived from Human NOS3 around the non-phosphorylation site of S1179.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	NOS3; Nitric oxide synthase; endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	NOS3


## Image




Western Blot analysis of A549 cells using NOS3 Polyclonal Antibody



Western Blot analysis of A549 cells using NOS3 Polyclonal Antibody



Western blot analysis of mouse-brain lysis using NOS3 antibody.



Western blot analysis of mouse-brain lysis using NOS3 antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.