




NTRK1 Antibody

Product Code	CSB-PA020170
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04629
Immunogen	Synthesized peptide derived from Human Trk A around the non-phosphorylation site of Y701.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NTRK1; MTC; TRK; TRKA; High affinity nerve growth factor receptor; Neurotrophic tyrosine kinase receptor type 1; TRK1-transforming tyrosine kinase protein; Tropomyosin-related kinase A; Tyrosine kinase receptor; Tyrosine kinase receptor A;
Immunogen Species	Homo sapiens (Human)
Target Names	NTRK1

Image



Western Blot analysis of 293T NIH-3T3 cells using Trk A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.